



THE SOLUTION FOR RAPID CONTROLLED DISINFECTION OF SMALL SURFACES AND EQUIPMENT



Klerwipe™ Impregnated Sterile Biocide Wipes

Ergonomic, simple to use pouch wipes with a wide range of applications

.....

These sterile 100% polyester or low-particulate pouch wipes are impregnated with a choice of Klercide Biocides, each provided in a convenient resealable pouch.

They are highly absorbent, allowing rapid disinfection of surfaces, and they remain wet throughout use, ensuring a consistent application of biocide to the whole surface area.

Klerwipe Biocide wipes:

- ▲ Are impregnated with a choice of efficacious biocides including sporicidal options
- ▲ Provide effective rotational use
- ▲ Feature Kler-tear™ easy open feature on outer packaging - ideal for gloved hands
- ▲ Have easy to use, resealable pouches ensuring the wipes don't dry out
- ▲ Are ideal for sterile environments



Klerwipe™ Impregnated Sterile Biocide Wipes



Ergonomic, simple to use pouch wipes with a wide range of applications

Why use Klerwipe Impregnated Sterile Biocide Wipes?

EFFECTIVE

All wipes are impregnated with 0.2 micron filtered fluids from the Klercide range of biocides.

They are double bagged for assured transfer and provide optimal effectiveness against a wide range of organisms commonly encountered in a cleanroom.

LOW RISK

The selected biocides are non-toxic and non-corrosive under normal use, so can be used in most cleanroom applications.

Being impregnated, Klerwipe pouch wipes eliminate the potential risk, mixing a concentrate or a dilution error and deliver the correct concentration of active ingredient to the surface being disinfected.

The Klerwipe pouch wipes are either packed in a Grade C (ISO Class 7) cleanroom and sterilised by gamma irradiation or aseptically processed, under a Grade A (ISO Class 5) laminar flow system to a validated procedure.

SIMPLE TO USE

The 200 x 200 mm wipes are made from a highly absorbent 100% polyester or non-woven, binder free, hydroentangled cellulose/polyester blend.

Each pouch is designed for sessional use and contains 10 or 15 wipes so there is no wastage at the end of a manufacturing session. The pouch is easily resealable to prevent the wipes drying out. The pouch itself does not generate any particulate contamination.

VALIDATION

Manufactured by an ISO 9001:2015 and ISO 14001:2015 certified company.

Certificates of Analysis and Sterility are available for every batch of Klerwipe pouch wipes. The shelf life of the product is dependent on the wipe selected.



Klerwipe Sterile Impregnated Biocide Pouch Wipes

WORLDWIDE HEADQUARTERS

1 Ecolab Place
St. Paul, MN 55102
USA

www.ecolab.com/lifesciences

EUROPE HEADQUARTERS

Richtistr. 7
8304 Wallisellen
Switzerland

www.ecolablifesciences.com

**Efficient.
Effective.
Convenient.**

APPLICATION AREA

For use on hard surfaces.

INSTRUCTIONS FOR USE

For optimum results ensure complete wetting of the surfaces. Re-seal when not in use. Dispose of packaging only if emptied and closed.

INGREDIENTS

See individual products for details.

PACKAGING INFORMATION

Klerwipe Impregnated Sterile Biocide Wipes

UNIT SIZE	UNIT OF SALE
200 x 200 mm pouch wipe	20 x 15 wipes*

*Klercide Low Residue Quat is 20 x 10 wipes

FOR FULL PRODUCT ORDER CODES, PLEASE VISIT OUR WEBSITE.

USE BIOCIDES SAFELY, ALWAYS READ THE LABEL AND PRODUCT INFORMATION BEFORE USE.

