

## Lab-Tips® - LTP70F

### Small Knitted Polyester Swab with Flexible Tip

Cleanroom laundered double layer 100% polyester knit swab

The LTP70F Swab is an excellent swab for cleaning grooved, slotted and recessed areas. Its cleanroom laundered knitted polyester head is extremely clean and durable. A compact handle and internal

flexible paddle provide maximum dexterity in confined areas.



#### Key Attributes

- Cleanroom laundered 100% continuous filament polyester knit head
- Short polypropylene handle with a flexible paddle head
- Thermally bonded head
- Laundered and packaged in Class 10 {ISO Class 4 } cleanroom

#### Benefits

- Critically low particles, fibers, ions and extractables
- High abrasion resistance
- Chemically compatible with IPA, acetone and other solvents
- Soft and non abrasive with no contaminating adhesives
- Good sorbency
- Autoclavable
- Free from Silicone, amides and DOP

#### Environmental

- Complies with RoHS and REACH requirements

#### Applications

- Designed for the highest level of contamination control in critical processing applications
- Cleaning of sensitive equipment
- Applying lubricants and other liquids
- Removing excess materials
- Cleaning with solvents

#### Physical Characteristics

<b>Head Material</b>	Knitted Polyester	<b>Handle Thickness</b>	2.2mm {0.086"}
<b>Head Width</b>	4mm {0.157"}	<b>Handle Length</b>	59mm {2.32"}
<b>Head Thickness</b>	2.5mm {0.098"}	<b>Total Swab length</b>	70mm {2.75"}
<b>Head Length</b>	11mm {0.433"}	<b>Head Bond</b>	Thermal
<b>Handle Material</b>	Polypropylene	<b>Handle Text</b>	Berkshire
<b>Handle Width</b>	2.2mm {0.086"}		

#### Other Knitted Polyester tip swabs

- LTP1465 Long handle
- LTP125 Large head
- LTP70R Rigid tip

[www.berkshire.eu.com](http://www.berkshire.eu.com)

Contact: Tel 44 (0) 870 757 2877  
enquiries@berkshire.uk.com

Americas	Tel 1 413 528 2602	Fax 1 413 528 2614	info@berkshire.com
Europe	Tel 44 (0) 870 757 2877	Fax 44 (0) 870 757 2878	enquiries@berkshire.uk.com
SE Asia	Tel 65 6252 4313	Fax 65 6252 4312	enquiries@berkshire.com.sg
Japan	Tel 81 3 5827 2380	Fax 81 3 5827 2382	master@berkshire.co.jp

**Technical Data:**

ATTRIBUTE		LTP70F	TEST METHOD
Non-Volatile Residue	IPA mg/swab	0.010	IEST-RP.CC004.3, Sec 7.1.2
Ions (µg/swab)	Chloride	0.058	IEST-RP.CC004.3, Sec 7.2.2
	Sulfate	0.16	
	Nitrate	ND*	
	Phosphate	ND*	
	Fluoride	ND*	
	Potassium	0.0015	
	Calcium	0.020	
	Sodium	0.0048	
	Magnesium	0.0046	

**Notes:**

- Data shown are typical values
- Due to differences in test methods applied and equipment utilized by different swab manufacturers, valid product comparisons may only be obtained through side-by-side testing in the same test facility, under similar conditions.
- \* ND\* = below the detection limit of the equipment.

**Order Information:**

Product	Part Number	Size	Swabs/Pk	Pks/Cs	Case size
Small Polyester w/ Flexible Tip	LTP70F.20	70mm	100	20	235 x 105 x 85mm


Other Berkshire products



Wipers



Mop Systems



Face Masks



Glove Liners



Documentation Systems



Swabs